




TD 8000

Cutlery polisher

Made in Germany



TECHNICAL DATA

dimensions (l x w x h):	64 x 79 x 75 cm
weight:	85 kg
installed load max.:	1050 W
power supply:	220-240 V – 50/60 Hz

YOUR ADVANTAGES AT A GLANCE

- ▶ **quick and efficient**
dries and polishes up to 8.000 cutlery pieces per hour in only one working step
- ▶ **clean and hygienic**
the cutlery looks like new and is sterilized by an UVC-lamp
- ▶ **cost-saving and economic**
low operating expenses and cost-savings due to the gain of time
- ▶ **noiseless and smooth running**
special spring system for smooth running
- ▶ **solid and durable**
high-grade materials and a solid construction for long-lasting quality
- ▶ **user-friendly and low-maintenance**
simple and safe handling due to the well-thought-out design
- ▶ **technical and modern**
innovative and safe technology
- ▶ **service-orientated and guaranteed**
we offer spare parts-guarantee and service



TRATAMIENTO DE RESIDUOS, S.L.
Polígono Ibaiondo, 16
Tel.: 943 331 424 - Fax: 943 331 549
Apartado 84
20120 HERNANI (Gipuzkoa) SPAIN
E-mail: str2000@str-2000.com
www.str-2000.com



The clean and wet cutlery is given directly from the dishwasher baskets into the TD 8000.



The cutlery is spirally transported upwards through a high-grade granular substance and is dried and polished. The installed ventilator prevents that granulate comes out of the machine.



The regular use of the TD 8000 keeps your cutlery shining.

The simple operation of TD 8000

Whether chromium, stainless steel or even high-grade silver – the clean and wet cutlery is given directly from the dishwasher into the cutlery polishing machine TD 8000.

Via the internal vibrating conveyor it is spirally transported upwards through a high-grade granular substance. This purely natural product dries and polishes your cutlery and keeps it shining– with just a single run. The granulate has to be exchanged every 4 – 8 weeks depending on the frequency the machine is used. The machine operates very quietly due to a special spring system.

Your cutlery will not only pass the inspection by your guests, it even meets the strict EU-Hygienic Guidelines: a special UVC-lamp sterilizes both, the cutlery and the granulate.

The constant use of the TD 8000 reduces the oxidation of cutlery in a remarkable way and makes it look noticeably better. As manufacturer we offer you a 10-years spare parts-guarantee and service.

Special technical features

- ▶ control of all components with digital display of error codes
- ▶ PC-interface for error read-out
- ▶ display showing when exchange of granulate or maintenance become due
- ▶ big double-spaced display, selectable from seven languages
- ▶ electronic temperature control for an optimisation of drying process
- ▶ function test of UVC-lamp



Cutlery polishing machine TD 8000



TRATAMIENTO DE RESIDUOS, S.L.

Polígono Ibaiondo, 16
 Tel.: 943 331 424 - Fax: 943 331 549
 Apartado 84
 20120 HERNANI (Gipuzkoa) SPAIN
 E-mail: str2000@str-2000.com
 www.str-2000.com